


Accessories for YMC Glass Columns



Accessories supplied with ECO and ECO^{PLUS} 10-15 mm ID

Connection Principle




Piston Thread: 1/4"-28G Ferrule Nut Tubing 1/16" (0.8 x 1.6)

Connection to Column 1/4"-28 G		K.P316
Connection to Column/ System 1/4"-28 G		JR-55050
Connection to System 10-32 HPLC (GE ÄKTA)		JR-5800
Connection to System M6 (GE ÄKTA)		K.P307

Accessories supplied with ECO and ECO^{PLUS} 25-50 mm ID

Connection Principle



Piston Thread: 1/4"-28G Ferrule Nut Tubing 1/8" (1.6 x 3.2)

Connection to Column 1/4"-28 G		K.P316
Connection to Column/ System 1/4"-28 G		JR-55051
Connection to System M6 (GE ÄKTA)		K.P307

Which Fitting Do I Have?

System

Fitting for the System
M6 (GE ÄKTA)



K.P207 1/16" Tubing
K.P307 1/8" Tubing

System

Fitting for the System
10-32 (GE ÄKTA)



JR-5800

Column

Lock Screw for the
Glass Column 1/4"-28 G



K.P316

Column/System

Fitting for the System/Glass Column
1/4"-28 G



JR-55050 1/16" Tubing
JR-55051 1/8" Tubing